



**KAT Implants, LLC**

T 603.427.0084  
E lab@katimplants.com  
15 Rye Street, Suite 115  
Portsmouth, NH 03801

LAB ONLY      Packaged by \_\_\_\_\_  
Date sent \_\_\_\_\_

Doctor \_\_\_\_\_  
Patient \_\_\_\_\_  
Address \_\_\_\_\_  
State \_\_\_\_\_ Zip \_\_\_\_\_  
Date Sent \_\_\_\_\_  
Due Date \_\_\_\_\_  
Allow 10 working days in-lab time from the time case is received + 2 days return shipping.


**Impression Coping TT Check List**  
 Use Adjustment Tool prior to each impression  
 Place coping on implant all the way down  
 Clean implant post  
 Use light and heavy body  
 Check if coping is secure inside impression  
 If coping moves, retake impression!

Specific Instructions:

SHADE \_\_\_\_\_  
TOOTH # (s) \_\_\_\_\_  
Abutment size(s) \_\_\_\_\_

Qty.  
 Emax  
 Zirconia Full Contour  
 Zirconia layered with Emax  
 Layered Pink Gingiva  
 Custom Abutment (milling is included)  
 Abutment Milling  
 Milled provisional crowns and bridges  
 Soft tissue (no charge if it is KAT Abutment and restoration is made in our lab)  
 Reduction coping (fabricated if prep adjustment is necessary due to insufficient occlusal reduction or undercuts. Will call to inform).

Torque Indexed Abutments to **15 Ncm**  
Torque Non-Indexed Abutments to **25 Ncm**



Additional Instructions

Prices and services are subject to change without prior notice.

All rush cases must be pre-scheduled by calling 603-427-0084 before the case is shipped. Additional fee applies. Shipping supplies and ground labels are available upon request. Expediated inbound shipping is available for an additional fee, please call for specific pricing.