



Implants			
Two-Piece Implants			
Item Name	Item Number	Item Image	Price
tom-KAT 2.5 x 10	10050		\$70
tom-KAT 2.5 x 12	10051		\$70
tom-KAT 2.5 x 14	10052		\$70
tom-KAT 3.1 x 10	10053		\$70
tom-KAT 3.1 x 12	10054		\$70
tom-KAT 3.1 x 14	10055		\$70
Implant 3.5 x 8	0013508		\$162
Implant 3.5 x 10	0013510		\$162
Implant 3.5 x 12	0013512		\$162
Implant 3.5 x 14	0013514		\$162
Implant 4.3 x 6	0014306		\$162
Implant 4.3 x 7	0014307		\$162
Implant 4.3 x 8	0014308		\$162
Implant 4.3 x 10	0014310		\$162
Implant 4.3 x 12	0014312		\$162
Implant 4.3 x 14	0014314		\$162
Implant 5.0 x 6	0015006		\$162
Implant 5.0 x 7	0015007		\$162
Implant 5.0 x 8	0015008		\$162
Implant 5.0 x 10	0015010		\$162
Implant 5.0 x 12	0015012	\$162	
Implant 6.5 x 7	0016507		\$162
Implant 6.5 x 8	0016508		\$162
One-Piece Implants			
One-Piece Implant 2.5 x 10	0012500		\$70
One-Piece Implant 2.5 x 12	0012502		\$70
One-Piece Implant 2.5 x 14	0012504		\$70
One-Piece Implant 3.0 x 10	0013000		\$70
One-Piece Implant 3.0 x 12	0013002		\$70
One-Piece Implant 3.0 x 14	0013004		\$70
bob-KAT 2.5 x 10	0012506		\$79
bob-KAT 2.5 x 12	0012507		\$79
bob-KAT 2.5 x 14	0012508		\$79
bob-KAT 3.1 x 10	0013006		\$79
bob-KAT 3.1 x 12	0013007	\$79	
bob-KAT 3.1 x 14	0013008	\$79	

Non-Indexed Abutments				
Contoured Abutments				
Item Name	Item Number	Item Image	Price	
Contoured Abutment 4.2 x 6.5	10010		\$76	
Contoured Abutment 4.2 x 7.2	10011		\$76	
Contoured Abutment 4.6 x 6.5	10013		\$76	
Contoured Abutment 4.6 x 7.2	10014		\$76	
Contoured Abutment 5.0 x 6.5	10016		\$76	
Contoured Abutment 5.0 x 7.2	10017		\$76	
Contoured Abutment 5.4 x 6.5	10019		\$76	
Contoured Abutment 5.4 x 7.2	10020		\$76	
Prepable Abutments				
Prepable Abutment 4.2 x 6.5	10031			\$76
Prepable Abutment 4.2 x 8.5	10032			\$76
Prepable Abutment 4.6 x 6.5	10033			\$76
Prepable Abutment 4.6 x 8.5	10034			\$76
Prepable Abutment 5.0 x 6.5	10035	\$76		
Prepable Abutment 5.0 x 8.5	10036	\$76		
Prepable Abutment 5.4 x 6.5	10037	\$76		
Prepable Abutment 5.4 x 8.5	10038	\$76		
Angled Abutments				
Angled Abutment 4.2 10°	002204			\$129
Angled Abutment 4.2 20°	002201			\$129
Angled Abutment 4.6 10°	002205			\$129
Angled Abutment 4.6 20°	002202			\$129
Angled Abutment 5.0 10°	002207		\$129	
Angled Abutment 5.0 20°	002208		\$129	
Healing Abutments				
Healing Abutment 4.2	0064246			\$43
Healing Abutment 4.7	0064746			\$43
Healing Abutment 5.1	0065146			\$43
Healing Abutment 5.5	0065546			\$43
Healing Abutment L/B	0064244	\$43		
Healing Abutment Plug (10 pack)	10131	\$11		
Cast-to Abutment IL				
Cast-to Abutment IL	003100			Call
Burn-out Abutment IL				
Burn-out Abutment IL	003102			\$48

Indexed Abutments				
Temporary Abutments, Crown				
Item Name	Item Number	Item Image	Price	
Temporary Abutment, Crown 4.2	002035		\$107	
Temporary Abutment, Crown 4.7	002036		\$107	
Temporary Abutments, Bridge				
Temporary Abutment, Bridge 4.2	10065		\$107	
Temporary Abutment, Bridge 4.7	10066		\$107	
Temporary Abutment, Bridge 5.1	10067		\$107	
Temporary Abutment, Bridge 5.5	10068		\$107	
Temporary Abutment, Bridge 5.5	10068		\$107	
CI Abutments				
CI Abutment 4.2	10101		\$105	
CI Abutment 4.6	10102		\$105	
CI Abutment 5.0	10103		\$105	
CI Abutment 5.4	10104		\$105	
CI Abutment 5.4	10104		\$105	
CI Abutment Cement-on Sleeves				
Sleeve Straight 4.2	10125		\$13	
Sleeve Straight 4.6	10126		\$13	
Sleeve Straight 5.0	10127		\$13	
Sleeve Straight 5.4	10128		\$13	
Sleeve Angled 4.2/4.6	10129		\$13	
Sleeve Angled 5.0/5.4	10130		\$13	
Denture Abutments and Housing				
ID™ Abutment				
Immediate Denture Abutment	003003			\$49
Locator Abutment				
Locator Abutment 2.5	01714		\$115	
Locator Abutment 3.5	01715		\$115	
Locator Abutment 4.5	01716		\$115	
Locator Processing Package (2-pack)	008002			\$59
Locator Processing Package (10-pack)	008003			\$266
Ball Abutments				
Ball Abutment 2	002104			\$43
Ball Abutment 3	002105			\$43
Ball Abutment 4	002106			\$43
Ball Abutment 5	002107			\$43
Ball Abutment 10°	002108	\$76		
Ball Abutment Housing				
Housing 4.2 with O-ring - used with bob-KAT implants	10080			\$25
Housing 4.7 with O-ring - used with Ball Abutments	002110		\$25	
O-ring 4.2 (10 pack)	10081		\$8	
O-ring 4.7 (10 pack)	002109		\$8	

# Instrumentation

## Instrumentation Highlights

- Use 1.5 and 2.0 mm drills at 1000 rpm with copious irrigation.
- First laser mark on 1.5 and 2.0 drills is at 8 mm.
- Use Final Drills at 100 rpm without irrigation.
- First laser mark on Final Drills is at 6 mm.
- Use Reamers at 100 rpm without irrigation.
- Reamer and Healing Abutment diameters have to be the same.
- Guide Pin Laser marks start at 6 mm; marks are placed every 2 mm.
- Implant Driver M is used with a Ratchet Wrench.
- Implant Driver M is designed to withstand forces up to 100 Ncm.
- Implant Driver HP is used with a Handpiece.
- Torque higher than 25 Ncm may result in Implant Driver HP fracture.
- Placement Instrument is used for attachment of all abutments.
- Removal Instrument is used to remove all abutments.
- Depth Gage laser marks start at 6 mm; marks are placed every 2 mm.
- Milling Analog should be used whenever extensive abutment preparation is required.
- Sonicate milled abutments in an alcohol bath for 3–5 minutes after each milling session.
- Torque Wrench should be treated with care.
- Use Torque Wrench to attach the abutments and to verify implant stability at time of insertion.
- Implant and Milling Analogs are compatible with all abutments of KAT System.

Item Name	Item Number	Item Image	Price
Twist Drill 1.5 mm	004115		\$35
Twist Drill 1.5 mm Short	004116		\$35
Twist Drill 2.0 mm	004120		\$35
Twist Drill 2.0 mm Short	004121		\$35
Final Drill 2.7 mm	004150		\$70
Final Drill 2.7 mm Short	004161		\$70
Final Drill 2.9 mm	004151		\$70
Final Drill 2.9 mm Short	004162		\$70
Final Drill 3.2 mm	004152		\$70
Final Drill 3.2 mm Short	004163		\$70
Final Drill 3.5 mm	004153		\$70
Final Drill 3.5 mm Short	004164		\$70
Final Drill 3.9 mm	004154		\$70
Final Drill 3.9 mm Short	004165		\$70
Final Drill 4.2 mm	004155		\$70
Final Drill 4.2 mm Short	004166		\$70
Final Drill 4.7 mm	004156		\$70
Final Drill 4.7 mm Short	004167		\$70
Final Drill 4.85 mm	004181		\$70
Final Drill 4.85 mm Short	004182		\$70
Final Drill 5.6 mm	004157		\$70
Final Drill 5.6 mm Short	004168		\$70
Final Drill 6.0 mm	004158		\$70
Final Drill 6.0 mm Short	004169		\$70
Trial Implants 2.5 - 6.5 mm	Specify Size		\$42
Drill Extender	004180		\$70
Guide Pin	004300		\$16
Retaining Screw	004009		\$16
Open Tray Impression Coping, 7mm	10120		\$16
Open Tray Impression Coping, 12mm	10121		\$16
Impression Coping TT 3.5	10117		\$32
Impression Coping TT 5.5	10118		\$32
Impression Coping TT 8.5	10119		\$32
Impression Coping TT Adjustment Tool	10111		\$27
Reamer 4.2 mm	004043		\$87
Reamer 4.7 mm	004048		\$87
Reamer 5.1 mm	004052		\$87
Reamer 5.5 mm	004056		\$87
Reamer L/B	004067		\$87
Square Driver FG	004011		\$76
Square Driver, Short	004010S		\$76
Square Driver, Long	004010L		\$76
Implant Driver M, Short	004003S		\$135
Implant Driver M, Long	004003L		\$135
Implant Driver HP, Short	004201		\$135
Implant Driver HP, Long	004203		\$135
Placement Instrument, Short	004001S		\$43
Placement Instrument, Long	004001L		\$43
Removal Instrument, Short	004002S		\$43
Removal Instrument, Long	004002L		\$43
Depth Gage	004012		\$76
Removal Screw, Short	004008S		\$43
Removal Screw, Long	004008L		\$43
Placement Screw, Short	004007S		\$43
Placement Screw, Long	004007L		\$43
Milling Analog	005003		\$87
Torque Wrench	009005		\$440
Instrumentation O-ring (10)	009100		\$11
Implant Analog	005001		\$38
Ball Abutment Placement Inst.	004021		\$76
Locator Placement Inst.	008015		\$76
Abutment Remover	0021059		\$134
Ball Abutment Removal Tip	0021058		\$81
Locator Removal Tip	008022		\$81
Casting Clean-up Drill	003105		\$21



### Aseptico Motors and Handpieces

Call for price.



### Surgical Kit

Includes all instrumentation to place and restore implants.

#009002  
\$2,500



### Restorative Kit

Includes:  
Placement/ Removal inst. S, L  
Placement/ Removal screw S, L  
Square Driver, S, L

#009003  
\$397



### Reamer Kit

Includes:  
Reamers 4.2, 4.7, 5.1, 5.5, L/B  
Placement/ Removal instrument S, L  
Square Driver, FG  
Impression TT adjustment tool

#009011  
\$552



### Laboratory Kit

Includes:  
Placement/ Removal inst. L  
Placement/ Removal screw L  
Square Driver, L

#009007  
\$161



### Short Drill Kit

Includes short length drills in  
1.5 - 6.0mm

#009009  
\$657